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Helene Derksen

### Literatur:

- Čelakovská J, Čermáková E, Andrýs C et al.: Sensitization to latex and food allergens in atopic dermatitis patients according to ALEX2 Allergy Explorer test. *Molecular Immunology* 2024; 175: 89-102  
<https://doi.org/10.1016/j.molimm.2024.09.002>
- Barrale M, Mazzucco W, Fruscione S et al.: An Assessment of the Currently Available Molecular Assay for the Diagnosis of Anisakis Sensitization. *International Journal of Molecular Sciences* 2025; 26 (7), 3033.  
<https://doi.org/10.3390/ijms26073033>
- Ukleja-Sokołowska N, Lis K, Bartuzi M et al.: Extended IgE profile of shrimp-sensitized patients based on multiplex examination ALEX2 – Allergy Explorer. *Advances in Dermatology and Allergology* 2023; 40 (5): 661-669.  
<https://oi.org/10.5114/ada.2023.132071>
- Kazancıoğlu A, Ocak M, Sahiner UM et al.: Natural history of sesame allergy in pediatric patients: Insight from a retrospective analysis. *Pediatric Allergy and Immunology* 2025; 36 (1): e70022.  
<https://doi.org/10.1111/pai.70022>
- Mederos-Luis E, Galán T, Martínez MJ et al.: Clinical and Immunological Features in Limpet (*Patella* sp.) Allergy in Subtropical Areas: A New Trigger for Food Anaphylaxis. *Allergies* 2025; 5 (1), 2.  
<https://doi.org/10.3390/allergies5010002>
- Hamilton RG, Holbreich M, Bronzert C et al.: Screening asthmatics for atopic status using the Allergy Explorer (ALEX2) macroarray. *Journal of Asthma* 2024; 61 (9): 1050-1057.  
<https://doi.org/10.1080/02770903.2024.2324839>

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